


# ELEVATED PLUS MAZE



Drawings and images are not shown proportionally for illustration purposes.

## Features & Benefits

- » Made of durable and easy to clean ABS plastic
- » Available in black, white or beige
- » One surface of the maze is smooth; the other is slightly textured
- » Designed for optimal use with SDI ANY-maze video tracking system that features:
  - › User defined zones
  - › Tracks time spent and distance traveled, speed and resting time in each zone
  - › Customizable Reports

## ELEVATED PLUS MAZE SPECIFICATIONS

	MAZE DIMENSIONS					
Mouse	(A) 26"	(B) 12"	(C) 2"	(D) 15.25"	(E) 21.25"	(F) .125" (lip)
Rat	(A) 43.5"	(B) 19.75"	(C) 4"	(D) 19.5"	(E) 31.5"	(F) .35" (lip)
Material Composition	.25" ABS plastic					
Color Options	Black, White, Beige					

## PRODUCT NAME

## PART NUMBERS

Plus Maze for Mice	Black 7001-0225	White 7001-0336	Beige 7001-0316
Plus Maze for Rats	Black 7001-0227	White 7001-0338	Beige 7001-0337

**Note:** Beige is the recommended color as it provides a neutral background for tracking both light and dark animals with the SDI ANY-maze Video Tracking System.