


RADIAL ARM MAZE



Drawings and images are not shown proportionally for illustration purposes.

Features & Benefits

- » Made of durable and easy to clean ABS plastic or acrylic
- » Available in black, white, beige or clear walls with colored or clear floors
- » Interchangeable plain and pellet feeder floors
- » Includes manual doors at each arm
- » Designed for optimal use with SDI ANY-maze video tracking system that features:
 - › User defined zones
 - › Tracks distance traveled, speed and resting time in each zone
 - › Customizable Reports

RADIAL ARM MAZE SPECIFICATIONS

	MAZE DIMENSIONS					
Mouse	(A) 26.45"	(B) 9"	(C) 8.45"	(D) 3"	(E) 6"	(H) .75" Diam. x .5" Deep
Rat	(A) 47.86"	(B) 18.5"	(C) 10.86"	(D) 4"	(E) 10"	(H) 1" Diam. x .75" Deep
Material Composition	.25" ABS Plastic - solid color mazes, .25" acrylic for clear wall mazes					
Color Options	Black, White, Beige or Clear Walls					

PRODUCT NAME

PART NUMBERS

Y Maze for Mice	Black 7001-0440	White 7001-0437	Beige 7001-0443	
Y Maze for Rats	Black 7001-0449	White 7001-0446	Beige 7001-0452	
Clear Walls for Mice	White Floors 7001-0438	Black Floors 7001-0441	Beige Floors 7001-0442	Clear Floors 7001-0444
Clear Walls for Rats	White Floors 7001-0447	Black Floors 7001-0450	Beige Floors 7001-0451	Clear Floors 7001-0453

Note: Beige is the recommended color as it provides a neutral background for tracking both light and dark animals with the SDI ANY-maze Video Tracking System.