


T MAZE



Drawings and images are not shown proportionally for illustration purposes.

Features & Benefits

- » Made of durable and easy to clean ABS plastic or acrylic
- » Available in black, white, beige or clear walls with colored or clear floors
- » Interchangeable plain and pellet feeder floors
- » Includes manual doors at each upper arm and starter door in base leg
- » Designed for optimal use with SDI ANY-maze video tracking system that features:
 - › User defined zones
 - › Tracks distance traveled, speed and resting time in each zone
 - › Customizable Reports

T MAZE SPECIFICATIONS

	MAZE DIMENSIONS						
Mouse	(A) 24"	(B) 3"	(C) 10.5"	(D) 15"	(E) 5"	(F) 3.5"	(H) .75" Diam. X .5" Deep
Rat	(A) 36"	(B) 4"	(C) 16"	(D) 22"	(E) 10"	(F) 7.25"	(H) 1" Diam. X .75" Deep
Material Composition	.25" ABS Plastic - solid color mazes, .25" acrylic for clear wall mazes						
Color Options	Black, White, Beige or Clear Walls						

PRODUCT NAME

PART NUMBERS

T Maze for Mice	Black 7001-0404	White 7001-0401	Beige 7001-0407	
T Maze for Rats	Black 7001-0413	White 7001-0410	Beige 7001-0416	
Clear Walls for Mice	White Floors 7001-0402	Black Floors 7001-0405	Beige Floors 7001-0406	Clear Floors 7001-0408
Clear Walls for Rats	White Floors 7001-0411	Black Floors 7001-0414	Beige Floors 7001-0415	Clear Floors 7001-0417

Note: Beige is the recommended color as it provides a neutral background for tracking both light and dark animals with the SDI ANY-maze Video Tracking System.