


# WATER TEST COMPONENTS

## Features & Benefits

- » Tank available in blue or white
- » Platforms are clear acrylic either fixed height or adjustable height
- » Walls are seamless, eliminating inadvertent cues
- » Designed for optimal use with SDI ANY-maze video tracking system that features:
  - › User defined zones
  - › Tracks distance traveled, speed and resting time in each zone
  - › Customizable Reports



Note: Drawings and images are not shown proportionally for illustration purposes.

## WATER TEST COMPONENTS SPECIFICATIONS

Mouse Tank Dimensions	(A) 48"	(B) 32"	235 GAL, 49 LBS (empty)				
Rat Tank Dimensions	(A) 72"	(B) 32"	525 GAL, 84 LBS, 73 LBS (empty)				
Platform Dimensions							
Square Fixed	(C) 12" x 12"	(D) 12.76"	(E) 4" x 4"	Round Fixed	(C) 12" x 12"	(D) 12.76"	(F) 4" diam.
Square Adjustable	(C) 12" x 12"	(D) 12.76"- 22"	(E) 4" x 4"	Round Adjustable	(C) 12" x 12"	(D) 12.76"- 22"	(F) 4" diam.
Platform Material Composition	Acrylic						
Tank Material Composition	Polyethylene						
Tank Color Options	Blue, White						

PRODUCT NAME	PART NUMBER	
Water Tank for Mice (4')	Blue 7000-0720	White 7000-0723
Water Tank for Rat (6')	Blue 7000-0718	White 7000-0719
Water Maze Platforms- Fixed Height	Clear Square Top 7500-0271	Clear Round Top 7500-0274
Water Maze Platforms- Adjustable Height	Clear Square Top 7500-0309	Clear Round Top 7500-0310
Drain Kit	7500-0217	