


Mazes & Enclosures

# Y MAZE



Drawings and images are not shown proportionally for illustration purposes.

## Features & Benefits

- » Made of durable and easy to clean ABS plastic or acrylic
- » Available in black, white, beige or clear walls with colored or clear floors
- » Interchangeable plain and pellet feeder floors
- » Includes manual doors at each upper arm and starter door in base leg
- » Designed for optimal use with SDI ANY-maze video tracking system that features:
  - › User defined zones
  - › Tracks distance traveled, speed and resting time in each zone
  - › Customizable Reports

## Y MAZE SPECIFICATIONS

	MAZE DIMENSIONS					
Mouse	(A) 3"	(B) 15"	(C) 15"	(D) 5"	(E) 3.5"	(H) .75" Diam. X .5" Deep
Rat	(A) 4"	(B) 22"	(C) 22"	(D) 10"	(E) 7"	(H) 1" Diam. X .75" Deep
Material Composition	.25" ABS Plastic - solid color mazes, .25" acrylic for clear wall mazes					
Color Options	Black, White, Beige or Clear Walls					

## PRODUCT NAME

## PART NUMBERS

Y Maze for Mice	Black 7001-0422	White 7001-0419	Beige 7001-0425	
Y Maze for Rats	Black 7001-0431	White 7001-0428	Beige 7001-0434	
Clear Walls for Mice	White Floors 7001-0420	Black Floors 7001-0423	Beige Floors 7001-0424	Clear Floors 7001-0426
Clear Walls for Rats	White Floors 7001-0429	Black Floors 7001-0432	Beige Floors 7001-0433	Clear Floors 7001-0435

**Note:** Beige is the recommended color as it provides a neutral background for tracking both light and dark animals with the SDI ANY-maze Video Tracking System.