


# ZERO MAZE



*Drawings and images are not shown proportionally for illustration purposes.*

## Features & Benefits

- » Made of durable and easy to clean ABS plastic
- » Available in black, white or beige
- » One surface of the maze is smooth; the other is slightly textured
- » Designed for optimal use with SDI ANY-maze video tracking system that features:
  - › User defined zones
  - › Tracks time spent and distance traveled, speed and resting time in each zone
  - › Customizable Reports

## ZERO MAZE SPECIFICATIONS

### MAZE DIMENSIONS

Mouse	Inside Wall: (A) 20" Diam. / 62.83" Circum.	Outside Wall: (B) 24" Diam. / 75.40" Circum. (C) 6"
	3 Legs Adjustable for leveling (not shown): 19.75" - 20.75"	
Rat	Inside Wall: (A) 40" Diam. / 125.66" Circum.	Outside Wall: (B) 48" Diam. / 150.8" Circum. (C) 8"
	6 Legs Adjustable for leveling (not shown): 19.75" - 20.75"	
Material Composition	.25" ABS Plastic	
Color Options	Black, White, Beige	

## PRODUCT NAME

## PART NUMBERS

Zero Maze for Mice	Black 2325-0230	White 2325-0229	Beige 2325-0231
Zero Maze for Rats	Black 2325-0227	White 2325-0226	Beige 2325-0228

**Note:** Beige is the recommended color as it provides a neutral background for tracking both light and dark animals with the SDI ANY-maze Video Tracking System.