

iWorx EMG Recording System

Record Four Channels of Surface EMG

The RS-EMG4 is a 4 Channel Biopotential Amplifier that allows multiple channels of Surface EMG to be recorded simultaneously from a single exercising or resting subject. The recorder is furnished with LabScribe™ Recording and Analysis Software and connects directly to a USB port on either a Macintosh or Windows PC.



Measure and Analyze

- Surface EMG
- Maximal voluntary contraction (MVC)
- Muscle force
- Fatigue
- Strength and Performance Testing

LabScribe Software

- Record four channels of EMG simultaneously
- Easy export
- Built in Journal for creating reports
- Lab Site license included
- Free software updates
- EMG analysis included

General

- Leads Color coded
- Enclosure Plastic
- Certification CE
- Power USB
- Warranty Protected by a 3 year warranty
- Software LabScribe™



iWorx Systems, Inc.