

iWorx EMG Recording System

Record Eight Channels of Surface EMG

The RS-EMG8 is an 8 Channel Biopotential Amplifier that allows multiple channels of Surface EMG to be recorded simultaneously from a single exercising or resting subject. The recorder is furnished with LabScribe™ Recording and Analysis Software and connects directly to a USB port on either a Macintosh or Windows PC.



Measure and Analyze

- Surface EMG
- Maximal voluntary contraction (MVC)
- Muscle force
- Fatigue
- Strength and Performance Testing

LabScribe Software

- Record eight channels of EMG simultaneously
- Easy export
- Built in Journal for creating reports
- Lab Site license included
- Free software updates
- EMG analysis included

General

- LeadsColor coded
- EnclosurePlastic
- CertificationCE
- PowerUSB
- WarrantyProtected by a 3 year warranty
- SoftwareLabScribe™



iWorx Systems, Inc.